



**MGL-FN-556/1~800mW**




**LD PUMPED ALL-SOLID-STATE  
YELLOW GREEN LASER  
AT 556nm**

Yellow green laser at 556nm is made features of ultra compact, long lifetime, low cost and high stability, which is used in scientific experiment, optical sensor, measurement, instrument, communication, spectrum analysis, etc.



**SPECIFICATIONS**

Wavelength (nm)	556±1	
Operating mode	CW	
Output power (mW)	>1, 5, 10, 20, ... ,500	>500, ... ,800
Power stability (rms, over 4 hours)	<1%, <2%, <3%	
Transverse mode	TEM <sub>00</sub>	Near TEM <sub>00</sub>
M <sup>2</sup> factor	<1.2	
Beam diameter at the aperture (1/e <sup>2</sup> , mm)	~2.0	
Beam divergence, full angle (mrad)	<1.5	
Polarization ratio	>100:1 Vertical±5 degree (Horizontal Optional)	
Warm-up time (minutes)	<10	
Pointing stability after warm-up (mrad)	<0.05	
Beam height from base plate (mm)	27.4	
Operating temperature (°C)	10~35	
Power supply (90-264VAC)	PSU-H-LED/PSU-H-FDA/PSU-SR	
Modulation option	TTL on/off, 1Hz-1KHz, 1KHz-10KHz, 10KHz-30KHz; and Analog modulation option	
Expected lifetime (hours)	10000	
Warranty	1 year	



**Note: The laser head needs to be used on a heat sink with good heat dissipation.**

MGL-FN-556	PSU-H-LED	PSU-H-FDA	PSU-SR
<p style="text-align: center;">197(L)×70(W)×50(H) mm<sup>3</sup>, 1.5 kg</p>	<p style="text-align: center;">277 (L) ×145(W) ×106 (H) mm<sup>3</sup>, 2.6 kg</p>	<p style="text-align: center;">275(L) ×145(W) ×104(H) mm<sup>3</sup>, 2.3 kg</p>	<p style="text-align: center;">188(L) ×145(W) ×83(H) mm<sup>3</sup>, 1.2kg</p>